



DDS CALORIMETERS

Scientific Analytical Calorimeter Solutions



CAL3K-A CALORIMETER

Oxygen Bomb Calorimeter

MANUFACTURING SUPERB CALORIMETERS FOR TODAY'S ANALYTICAL NEEDS

www.ddscalorimeters.com

CAL3K-A CALORIMETER

The CAL3K-A is the second in the range of innovative new oxygen bomb calorimeters from DDS. It is the companion product to the CAL3K-AP, but with external oxygen filling. The new range, from the engineers who designed the CAL2K oxygen bomb calorimeter range, features higher speeds and improved accuracy with a small carbon footprint (uses less resources; no water required, lower power consumption). The new bayonet bomb vessel with the option of the adiabatic, isothermal or dynamic calorimetry methods results in optimized speed and accuracy. This can be tailored to your laboratory needs by parameter settings via the USB port. The CAL3K-A conforms to ASTM D5865-12, ASTM D4809-13, ISO 1928-2009, and DIN 51900-2 international standards. The traditional CAL2K-4 bomb vessel has been replaced with a self-locking and sealing bayonet vessel with built-in temperature sensors. The vessel is completely removable for easy sample management, cooling, cleaning and maintenance.

COMPLETE SYSTEM

The oxygen bomb calorimeter, filling station and air cooler are operated together for effective routine sample determination, using 2 or more bomb vessels. It is recommended to use the air cooler for optimum results and faster throughput.



COMPLETE SYSTEM

FOR USE WITH THE CAL3K-A.

The CAL3K-A is best suited for a production environment with shared oxygen supply (to other users) and robust handling.

The following accessories can be added to the calorimeter:

- CAL3K Air Cooler
- CAL3K Manual Bayonet Bomb Vessel
- **Analytical Balance (sold separately)**
- High Pressure Oxygen Regulator

CAL³K
NEXT GENERATION CALORIMETERS

The complete CAL3K-A oxygen bomb calorimeter system contains all the parts and consumables necessary to set up the unit. The installation kits included with the setup of the calorimeter contains consumables for approximately 200 samples, depending on the type of sample being analyzed (coal analysis, animal feed analysis). Other samples like oil, might use more consumables as they are corrosive and could cause wear and tear. Additional consumables can be purchased separately from our authorized agents.

The complete calorimeter system is delivered with installation kits for the calorimeter, filling station, air cooler and vessels (for approx. 200 firings). The CAL3K-A uses a combination of the Isothermal and Adiabatic and Dynamic (Isoperibol) calorimetry methods, while still using the dry method, i.e. it's waterless. The sample repeat speed is between 4-5 minutes.

The vessel is manually filled with oxygen via the external oxygen bomb filling station.

An average of 10 or more tests per hour can be achieved with 2 vessels.

CAL3K-A CALORIMETER

ADVANCED CAL3K-A FEATURES



TEMPERATURE CONTROL

No temperature control of room/lab required



FAULT FINDING

Extensive fault finding and testing



STEP-BY-STEP HELP

Screen prompts assist with step-by-step instructions to operate the calorimeter



7 CALIBRATION FIELDS

For different mode and different calorimeters



TEMPERATURE ACCURACY

Temperature accuracy of 10ppm (parts per million)(0.00001°C)



BAYONET BOMB VESSEL

Self-Locking and Self-Sealing manual bayonet bomb vessel



HIGH SPEED DETERMINATIONS

Choose between faster or more accurate determinations



RESULTS

Results in MJ/Kg, BTU/lb or Cal/g



AIR COOLER

No water required to cool the bomb vessel



EXTREMELY ACCURATE

Extremely accurate (%RSD - 0.05%) determination eliminates multiple sample repeats



COMPENSATION*

Compensation for firing energy and sulphur



BALANCE INTERFACE

Balance interface with baud speed setting



MULTIPLE COMM CHANNELS

1 Wired (USB) and 2 optional wireless channels



PRESET FIELDS

One default setting per mode



IMPROVED INITIAL TIMING

Based on drift and time, or time after drift



USER FRIENDLY

User Friendly Operation



OPERATING PARAMETERS

Operating parameters can be changed via USB interface in experimental mode



FULL LIMS SUPPORT

For the assignment, scheduling, and tracking of samples.



10 CALIBRATION AVERAGE

For variable amount of calibration average to suit your application



RESTRICT ACCESS

Operating parameter access is password restricted



FILTER DATA ON EVENTS

Extensive testing and detailed data viewing



LOW POWER CONSUMPTION

Very low power consumption. No temperature controlling required.



LARGE STORAGE

Up to 1024 results storage



ECO FRIENDLY

Eco Friendly - small carbon footprint. No water, low power consumption.



INTELLIGENT VESSEL

Intelligent vessel with built-in temperature sensing



TEMPERATURE RANGE

Extensive temperature range from 0°C to 70°C.



LINEAR SENSORS

Linear temperature sensing with platinum sensors



EVENT LOGGING

Built-in event logging for ~6000 events



SAFETY

Safety checks guarantee the safety of the operator.



NO WATER REQUIRED

No Water Bucket. No Spillage. No Measuring.



THREE OPERATING MODES

Isothermal, Adiabatic and ISOBATIC (and experimental)



MANUAL OXYGEN FILLING

Makes use of an external oxygen filling station

CAL3K-A CALORIMETER

The CAL3K Bomb Calorimeter Installation 3K-A-Kit includes :

- Power Supply (External 12V/3A)(3K-1-062)
- Mains Cable with plug and kettle connector (1.8m)(3K-1-122)
- Balance Cable (3K-1-117)
- PC Keyboard PS2 (3K-1-061)
- USB Communication Cable (3K-1-084)
- Preparation Stand (3K-4-49)
- Tweezers (3K-1-081)
- Certified Benzoic Acid Tablets (150 x 0.5g tablets per bottle)(3K-4-71)
- Printed Installation Manual (INSTALLATION-MANUAL)
- Vessel Spanner (3K-4-72)
- USB 32Gb Green Memory Flash Drive (3K-1-043)

CAL3K-A KIT AND CONSUMABLES

The CAL3K Air Cooler Installation 3K-2-Kit Includes :


- Power Supply (External 12V/1.25A direct 2 PIN) (3K-1-005)

The CAL3K Bayonet Vessel Installation 3K-5-Kit Includes :

- Bottom Centre Electrode (3K-4-36)
- Top Centre Electrode (3K-4-46)
- Base Nozzle O-Ring (3K-4-23)
- Outside Electrode (3K-4-37)
- Single Crucibles (3K-4-47)
- Deflector Plate (3K-4-05)
- Firing Cotton (Bundle of 100 threads)(3K-4-65)
- Firing Wire (1 packet of 5 wires)(3K-4-32)
- Large Lid O-Ring (3K-4-16)
- Top and Bottom O-Ring in Vessel Lid (3K-4-22)


The CAL3K Filling Station Installation 3K-3-Kit Includes :

- Nozzle O-Ring for Defiller Cap and Filling Station (3K-3-29)
- Nozzle O-Ring Used in 3K-3 Jet Assembly (3K-1-080)
- Oxygen Regulator Connection Kit (3K-3-21)
- Defiller Cap (3K-3-22)
- Emergency Deflate Cap (3K-3-18)
- High Pressure Pipe 4mm (White)(PIPE-4MM-CLEAR)



CAL3K Bayonet Bomb Vessel & Air Cooler

TECHNICAL SPECIFICATIONS



LCD Display

Large display for easy viewing.

x3

Easy to Use

With 3 different operating modes



Compact Size

Approximately 350mm x 280mm x 240mm



Lightweight

Light weight for easy moving.
Approximately 12kg.



TUV CE Certification

Complies with ASTM, DIN and ISO International Standards.

Specification	Information
Working (Operating) Temperature	15-70°C
Storage Temperature	0-70°C
Temperature Resolution	0.000001°C
Reproducibility/Repeatability	<0.01% RSD
Resolution	0.0001 MJ/Kg
Results per hour	12 samples per hour using 3 bayonet bomb vessels
Measuring range max.	99MJ, 99000J
Working temperature min.	0°C
Working temperature max.	50°C
Temperature Measurement Resolution	10ppm (parts per million)
Cooling Medium	Air

TECHNICAL SPECIFICATIONS

Specification	Information
Type of Cooling	Airflow
Oxygen Operating Pressure Max	40 bar
Balance/Scale Interface	RS232
Printer Interface	RS232
PC Interface	USB
Power Input	12W
Interface External Keyboard	PS2
Oxygen Filling	Manual
De gasification	Manual
Halogen (Decomposition) Vessel	Yes, optional
Analysis according to DIN 51900	Yes
Analysis according to ASTM D240	Yes
Analysis according to ASTM D4809	Yes
Analysis according to ASTM D5865	Yes
Analysis according to ASTM E711	Yes
Analysis according to ISO 1928	Yes
Dimensions	350mm x 280mm x 240mm
Weight	~ 12.000kg
Permissible Ambient Temperature	0-35°C
Permissible Relative Humidity	80%
RS232 Interface	Yes
USB Interface	Yes
Voltage	220-240 / 100-120V, 12VDC, 3Amp
Frequency	50/60 Hz

Please Note : Technical Specifications subject to change without prior notice.

Please contact our team for accurate technical specifications at the time.





DDS CALORIMETERS

Scientific Analytical Calorimeter Solutions



CAL3K-AP CALORIMETER

Oxygen Bomb Calorimeter

MANUFACTURING SUPERB CALORIMETERS FOR TODAY'S ANALYTICAL NEEDS

www.ddscalorimeters.com

CAL3K-AP CALORIMETER

The CAL3K-AP is the first in the range of innovative new oxygen bomb calorimeters from DDS. The new range, from the engineers who designed the CAL2K oxygen bomb calorimeter range, features higher speeds and improved accuracy with a small carbon footprint (uses less resources; no water required, lower power consumption). The new bayonet bomb vessel with the option of the adiabatic, isothermal or dynamic calorimetry methods results in optimized speed and accuracy. This can be tailored to your laboratory needs by parameter settings via the USB port. The CAL3K-AP conforms to ASTM D5865-12, ASTM D4809-13, ISO 1928-2009, and DIN 51900-2 international standards. The traditional CAL2K-4 bomb vessel has been replaced with a self-locking and sealing bayonet vessel with built-in temperature sensors. The vessel is completely removable for easy sample management, cooling, cleaning and maintenance. The CAL3K can interface with multiple devices required for sample determination.

COMPLETE SYSTEM

The oxygen bomb calorimeter and air cooler are operated together for effective routine sample determination, using 2 or more bomb vessels. It is recommended to use the air cooler for optimum results and faster throughput.



COMPLETE SYSTEM

FOR USE WITH THE CAL3K-AP.

The CAL3K-AP is best suited for a clean and well maintained laboratory with a requirement for very fast and accurate results.

The following accessories can be added to the calorimeter:

- CAL3K Air Cooler
- CAL3K Auto Bayonet Bomb Vessel
- Analytical Balance (Sold Separately)

CAL³K
NEXT GENERATION CALORIMETERS

The complete CAL3K-AP oxygen bomb calorimeter system contains all the parts and consumables necessary to set up the unit. The installation kits included with the setup of the calorimeter contains consumables for approximately 200 samples, depending on the type of sample being analyzed (coal analysis, animal feed analysis). Other samples like oil, might use more consumables as they are corrosive and could cause wear and tear. Additional consumables can be purchased separately from our authorized agents.

The complete calorimeter system is delivered with installation kits for the calorimeter, air cooler and vessels (for approx. 200 firings). The CAL3K-AP uses a combination of the Isothermal and Adiabatic calorimetry methods, while still using the dry method, i.e. it's waterless. The sample repeat speed is between 4-6 minutes.

The vessel is automatically filled with oxygen inside the calorimeter, and defilled on completion of a determination.

An average of 8 or more tests per hour can be achieved with 2 vessels.

CAL3K-AP CALORIMETER

ADVANCED CAL3K-AP FEATURES



VARIABLE FILLING PRESSURE
Allows for variable filling pressure



FAULT FINDING
Extensive fault finding and testing



STEP-BY-STEP HELP
Screen prompts assist with step-by-step instructions to operate the calorimeter



7 CALIBRATION FIELDS
For different mode and different calorimeters



TEMPERATURE ACCURACY
Temperature accuracy of 10ppm (parts per million)(0.00001°C)



BAYONET BOMB VESSEL
Self-Locking and Self-Sealing bayonet bomb vessel



HIGH SPEED DETERMINATIONS
Choose between faster or more accurate determinations



RESULTS
Results in MJ/Kg, BTU/lb or Cal/g



AIR COOLER
No water required to cool the bomb vessel



EXTREMELY ACCURATE
Extremely accurate (%RSD - 0.05%) determination eliminates multiple sample repeats



COMPENSATION*
Compensation for firing energy and sulphur



BALANCE INTERFACE
Balance interface with baud speed setting



MULTIPLE COMM CHANNELS
1 Wired (USB) and 2 optional wireless channels



PRESET FIELDS
One default setting per mode



IMPROVED INITIAL TIMING
Based on drift and time, or time after drift



USER FRIENDLY
User Friendly Operation



OPERATING PARAMETERS
Operating parameters can be changed via USB interface in experimental mode



FULL LIMS SUPPORT
For the assignment, scheduling, and tracking of samples.



10 CALIBRATION AVERAGE
For variable amount of calibration average to suit your application



RESTRICT ACCESS
Operating parameter access is password restricted



FILTER DATA ON EVENTS
Extensive testing and detailed data viewing



LOW POWER CONSUMPTION
Very low power consumption. No temperature controlling required.



LARGE STORAGE
Up to 1024 results storage



TEMPERATURE CONTROL
No temperature control of room/lab required



ECO FRIENDLY
Eco Friendly - small carbon footprint. No water, low power consumption.



INTELLIGENT VESSEL
Intelligent vessel with built-in temperature sensing



LEAK DETECTION
Detects leaks and aborts the firing in the vent of a leak



TEMPERATURE RANGE
Extensive temperature range from 0°C to 70°C.



LINEAR SENSORS
Linear temperature sensing with platinum sensors



EVENT LOGGING
Built-in event logging for ~6000 events



SAFETY
Safety checks guarantee the safety of the operator.



NO WATER REQUIRED
No Water Bucket. No Spillage. No Measuring.



THREE OPERATING MODES
Isothermal, Adiabatic and ISOBATIC (and experimental)



AUTOMATIC OXYGEN FILLING
Automatic oxygen filling and de-filling of the reaction vessel

CAL3K-AP CALORIMETER

The CAL3K Bomb Calorimeter Installation 3K-AP-Kit includes :

- Power Supply (External 12V/3A)(3K-1-062)
- Mains Cable with plug and kettle connector (1.8m)(3K-1-122)
- Balance Cable (3K-1-117)
- USB Communication Cable (3K-1-084)
- PC Keyboard (PS2)(3K-1-061)
- Preparation Stand (3K-4-49)
- Tweezers (3K-1-081)
- Certified Benzoic Acid Tablets (150 x 0.5g tablets per bottle)(3K-4-71)
- PVC Flexible Tubing 6.4mm (3K-6-012)
- Exhaust O-Ring (3K-6-031)
- Oxygen Regulator Connection Kit (3K-3-21)
- High Pressure Pipe 4mm (White)(PIPE-4MM-CLEAR)
- Nozzle O-Ring using in 3K-3 Jet Assembly (3K-1-080)
- Filling Nut (3K-6-018)
- Inflation Nipple on Cap (3K-4-39)
- O-Ring Lubrication Grease Tube (3K-1-086)
- Printed Installation Manual (INSTALLATION-MANUAL)
- Fill Nut Tool (Complete)(3K-6-017)
- Fill Head O-Ring (3K-6-025)
- USB 32Gb Green Memory Flash Drive (3K-1-043)
- Pressure Leak Test Jig (3K-1-115)
- Compression Pipe fitting double ferrule 4mm (3K-6-033)


CAL3K-AP KIT AND CONSUMABLES

The CAL3K Air Cooler Installation 3K-2-Kit includes :

- Power Supply (External 12V/1.25A direct 2 PIN) (3K-1-005)


The CAL3K Bayonet Vessel Installation 3K-4-Kit includes :

- Nozzle O-Ring for Defiller Cap and Filling Station (3K-3-29)
- Outside Electrode (3K-4-37)
- Bottom Centre Electrode (3K-4-36)
- Single Crucibles (3K-4-47)
- Deflector Plate (3K-4-05)
- Firing Wire (1 packet of 5 wires)(3K-4-32)
- Firing Cotton (Bundle of 100 threads)(3K-4-65)
- Large Lid O-Ring (3K-4-16)



CAL3K Bayonet Bomb Vessel & Air Cooler

TECHNICAL SPECIFICATIONS




LCD Display

Large display for easy viewing.

x3

Easy to Use

With 3 different operating modes



Compact Size

Approximately 350mm x 280mm x 240mm



Lightweight

Light weight for easy moving.
Approximately 12kg.



TUV CE Certification

Complies with ASTM, DIN and ISO International Standards.

12 samples per hour using 3 bayonet bomb vessels

Specification	Information
Working (Operating) Temperature	15-70°C
Storage Temperature	0-70°C
Temperature Resolution	0.000001°C
Reproducibility/Repeatability	<0.01% RSD
Resolution	0.0001 MJ/Kg
Results per hour	9 samples per hour using 2 bayonet bomb vessels
Measuring range max.	99MJ, 99000J
Working temperature min.	0°C
Working temperature max.	50°C
Temperature Measurement Resolution	10ppm (parts per million)
Cooling Medium	Air

TECHNICAL SPECIFICATIONS

Specification	Information
Type of Cooling	Airflow
Oxygen Operating Pressure Max	40 bar
Balance/Scale Interface	RS232
Printer Interface	RS232
PC Interface	USB
Power Input	12W
Interface External Keyboard	PS2
Oxygen Filling	Automatic
De gasification	Automatic
Halogen (Decomposition) Vessel	Yes, optional
Analysis according to DIN 51900	Yes
Analysis according to ASTM D240	Yes
Analysis according to ASTM D4809	Yes
Analysis according to ASTM D5865	Yes
Analysis according to ASTM E711	Yes
Analysis according to ISO 1928	Yes
Dimensions	350mm x 280mm x 240mm
Weight	~ 12.000kg
Permissible Ambient Temperature	0-35°C
Permissible Relative Humidity	80%
RS232 Interface	Yes
USB Interface	Yes
Voltage	220-240 / 100-120V, 12VDC, 3Amp
Frequency	50/60 Hz

Please Note : Technical Specifications subject to change without prior notice.

Please contact our team for accurate technical specifications at the time.





DDS CALORIMETERS

Scientific Analytical Calorimeter Solutions



CAL3K-F CALORIMETER

Oxygen Bomb Calorimeter

MANUFACTURING SUPERB CALORIMETERS FOR TODAY'S ANALYTICAL NEEDS
www.ddscalorimeters.com

CAL3K-F CALORIMETER

The CAL3K-F is the fourth in the range of innovative new oxygen bomb calorimeters from DDS Calorimeters. The new range, from the engineers who designed the CAL2K Oxygen Bomb Calorimeter range, features higher speeds, improved accuracy and a small footprint with the now legendary resource reduction (no water required, lower power consumption, low maintenance).

The CAL3K-F Oxygen Bomb Calorimeter is highly expandable to cater for increasing laboratory work. Start off with the most economical single vessel system and expand it as the need arises to a system of 9 samples per hour. The system uses the proven bayonet vessel common to all CAL3K systems, and works in the dynamic (Isothermal) calorimetry mode. It has three communication ports with a host of outstanding features that caters for every application. The CAL3K interfaces with different devices, including PC, Balance, and Printer.

COMPLETE SYSTEM

The oxygen bomb calorimeter, filling station and air cooler are operated together for effective routine sample determination, using 2 or more bomb vessels. It is recommended to use the air cooler for optimum results and faster throughput.



COMPLETE SYSTEM

FOR USE WITH THE CAL3K-F.

The CAL3K-F is a flexible system for low to high throughput without compromising accuracy and repeatability. It is used in Food/Feed Analysis, Alternative Energy, Scientific Determination, Coal and Oil Production and Research and Quality Assurance. In short : wherever the calorific value of a solid/liquid sample must be determined.

The following accessories can be added to the calorimeter :

- Sartorius Balance (sold separately)
- RS232 Serial Printer (sold separately)
- Laptop / Stand-Alone PC (sold separately)
- High Pressure Oxygen Regulator (Requirement)

CAL³K
NEXT GENERATION CALORIMETERS

The complete CAL3K-F oxygen bomb calorimeter system contains all the parts and consumables necessary to set up the unit. The installation kits included with the setup of the calorimeter contains consumables for approximately 200 samples, depending on the type of sample being analyzed (coal analysis, animal feed analysis). Other samples like oil, might use more consumables as they are corrosive and could cause wear and tear. Additional consumables can be purchased separately from our authorized agents.

The complete calorimeter system is delivered with installation kits for the calorimeter, filling station, air cooler and vessels (for approx. 200 firings). The CAL3K-F uses Dynamic (Isothermal) calorimetry method, while still using the dry method, i.e. it's waterless. The sample repeat speed is between 5-10 minutes.

The vessel is manually filled with oxygen via the external oxygen bomb filling station.

An average of 2-9 tests per hour can be achieved depending on system option purchased.

CAL3K-F CALORIMETER

ADVANCED CAL3K-F FEATURES



TEMPERATURE CONTROL

No temperature control of room/lab required



FAULT FINDING

Extensive fault finding and testing



STEP-BY-STEP HELP

Screen prompts assist with step-by-step instructions to operate the calorimeter



7 CALIBRATION FIELDS

For different mode and different calorimeters



TEMPERATURE ACCURACY

Temperature accuracy of 10ppm (parts per million)(0.00001°C)



BAYONET BOMB VESSEL

Self-Locking and Self-Sealing manual bayonet bomb vessel



EXTREMELY ACCURATE

Extremely accurate (%RSD - 0.05%) determination eliminates multiple sample repeats



RESULTS

Results in MJ/Kg, BTU/lb or Cal/g



AIR COOLER

No water required to cool the bomb vessel



MULTIPLE COMM CHANNELS

1 Wired (USB) and 1 optional wireless channels



COMPENSATION*

Compensation for firing energy and sulphur



BALANCE INTERFACE

Balance interface with baud speed setting



USER FRIENDLY

User Friendly Operation



PRESET FIELDS

One default setting per mode



OPERATING PARAMETERS

Operating parameters can be changed via USB interface in experimental mode



FULL LIMS SUPPORT

For the assignment, scheduling, and tracking of samples.



10 CALIBRATION AVERAGE

For variable amount of calibration average to suit your application



RESTRICT ACCESS

Operating parameter access is password restricted



LARGE STORAGE

Up to 1024 results storage



ECO FRIENDLY

Eco Friendly - small carbon footprint. No water, low power consumption.



INTELLIGENT VESSEL

Intelligent vessel with built-in temperature sensing



EVENT LOGGING

Built-in event logging for ~6000 events



LINEAR SENSORS

Linear temperature sensing with platinum sensors



SAFETY

Safety checks guarantee the safety of the operator.



NO WATER REQUIRED

No Water Bucket. No Spillage. No Measuring.



MANUAL OXYGEN FILLING

Makes use of an external oxygen filling station

CAL3K-F CALORIMETER

The CAL3K-F Bomb Calorimeter Installation Kit includes :

- Power Supply 12V/1.25A (3K-1-055)
- Balance Cable (3K-1-117)
- PC Cable USB (3K-1-084)
- PC Keyboard PS2 (3K-1-061)
- Preparation Stand (3K-4-49)
- Stainless Steel Tweezers (3K-1-081)
- Certified Benzoic Acid Tablets (150 x 0.5g tablets per bottle)(3K-4-71)
- Vessel Spanner for 3K-5) (3K-4-72)
- Printed Installation Guide (INSTALLATION-MANUAL)
- USB Memory Stick (including manuals, software, application notes etc) (3K-1-043)

CAL3K-F KIT AND CONSUMABLES

The CAL3K Air Cooler Installation Kit includes :


- Power Supply (3K-1-055)

The CAL3K Bayonet Vessel Installation Kit includes :

- 2 x Center Electrode Bottom Only (3K-4-36)
- 1 x Centre Electrode Top Only (3K-4-46)
- 2 x Nozzle O-Ring (3K-4-23)
- 2 x Outside Electrode (3K-4-37)
- 2 x Stainless Steel Crucibles (3K-4-47)
- 2 x Deflector Plate (3K-4-05)
- 1 x Packet Firing Cotton (100 threads) 3K-4-65
- 1 x Packet Firing Wire (5 wires) 3K-4-32
- 5 x Large Lid O-ring (3K-4-16)
- 10 x Top & Bottom O-ring (3K-4-22)


The CAL3K Filling Station Installation Kit includes :

- 5 x Nozzle O-ring (3K-3-29)
- 2 x Nozzle Piston O-ring (3K-1-080)
- 1 x Oxygen Regulator Kit (3K-3-21)
- 1 x Defiller Cap (3K-3-22)
- 1 x Emergency Deflate Cap (3K-3-18)



CAL3K Bayonet Bomb Vessel & Air Cooler

TECHNICAL SPECIFICATIONS



LCD Display

Large display for easy viewing.

x3

Easy to Use

With 3 different operating modes



Compact Size

Approximately 350mm x 280mm x 240mm



Lightweight

Light weight for easy moving.
Approximately 12kg.



TUV CE Certification

Complies with ASTM, DIN and ISO International Standards.

Specification	Information
Working (Operating) Temperature	15-70°C
Storage Temperature	0-70°C
Temperature Resolution	0.000001°C
Reproducibility/Repeatability	<0.01% RSD
Resolution	0.0001 MJ/Kg
Results per hour	Up to 8-9 samples per hour using 3 bayonet bomb vessels and 2 Air Coolers
Measuring range max.	99MJ, 99000J
Working temperature min.	1°C
Working temperature max.	50°C
Temperature Measurement Resolution	10ppm (parts per million)
Cooling Medium	Air

TECHNICAL SPECIFICATIONS

Specification	Information
Type of Cooling	Airflow
Oxygen Operating Pressure Max	40 bar
Balance/Scale Interface	RS232, 1200 to 38400 Baud (settable)
Printer Interface	RS232, 1.2Kb to 115.2Kb
PC Interface	USB
Power Input	2.4W
Interface External Keyboard	PS2
Oxygen Filling	Manual
De gasification	Manual
Halogen (Decomposition) Vessel	Yes, optional
Analysis according to DIN 51900	Yes
Analysis according to ASTM D240	Yes
Analysis according to ASTM D4809	Yes
Analysis according to ASTM D5865	Yes
Analysis according to ASTM E711	Yes
Analysis according to ISO 1928	Yes
Dimensions	350mm x 280mm x 240mm
Weight	~ 12.000kg
Permissible Ambient Temperature	1-35°C
Permissible Relative Humidity	80%
RS232 Interface	Yes
USB Interface	Yes
Voltage	220-240 / 100-120V, 12VDC, 1Amp
Frequency	50/60 Hz

Please Note : Technical Specifications subject to change without prior notice.


Please contact our team for accurate technical specifications at the time.





DDS CALORIMETERS

Scientific Analytical Calorimeter Solutions



CAL3K-S CALORIMETER

Oxygen Bomb Calorimeter

MANUFACTURING SUPERB CALORIMETERS FOR TODAY'S ANALYTICAL NEEDS

www.ddscalorimeters.com

CAL3K-S CALORIMETER

CAL3K-S is the fifth in the range of innovative new oxygen bomb calorimeters from DDS Calorimeters. The new range, from the engineers who designed the CAL2K Oxygen Bomb Calorimeter range, is a totally dry type calorimeter, making it environmentally friendly and energy efficient. The CAL3K-S also boasts the smallest carbon footprint in its class.

The CAL3K-S is an economical system and is typically used in applications where low sample analysis is required, like in Food/Feed Analysis, Alternative Energy, Manufacturing, Coal and Oil Production, Research, Universities, and Quality Assurance. In short: wherever the calorific value of a solid or liquid sample must be determined, the CAL3K-S is the economical calorimeter of choice.

COMPLETE SYSTEM

The CAL3K-S system includes a calorimeter, 1 x standard thread bomb vessel, and oxygen filling station. A spare bomb vessel (3K-5S) can be purchased, however this will not increase the speed of the calorimeter system.

The calorimeter system can complete 1 sample per hour. Part of the CAL3K-S system is the external oxygen filling station (3K-3), standard across all manual filling calorimeter systems.



COMPLETE SYSTEM

FOR USE WITH THE CAL3K-S.

The CAL3K-S is your entry level system suitable for low throughput without compromising on accuracy and repeatability. It is ideal for use in Food and Nutrition, Animal Feed Analysis, Universities and Research Departments, Alternative Energy Research, Waste and Waste Product Analysis.

The following are included:

- 3K-3 Manual Oxygen Filling Station
- 3K-5S Standard Bomb Vessel

The following can be optionally added to the system:

- Analytical Balance (sold separately)
- High Pressure Oxygen Regulator

The complete CAL3K-S oxygen bomb calorimeter system contains all the parts and consumables necessary to set up the system. The installation kits included with the setup of the calorimeter contain consumables for approximately 200 samples, depending on the type of sample being analysed, such as coal, food, or feed samples.

Other samples, such as oil, might use more consumables as they are corrosive and could cause wear and tear.

Additional consumables can be purchased separately from DDS or an authorised agent.

The sample repeat speed is 30 minutes. With Built-In Fan Cooler.

The vessel is manually filled with oxygen via the external manual oxygen filling station (3K-3).

1 Vessel system. Not expandable.

The system can analyse 1 sample per hour; a total of 8 samples per day.

CAL³K
NEXT GENERATION CALORIMETERS

CAL3K-S CALORIMETER

ADVANCED CAL3K-S FEATURES



TEMPERATURE CONTROL

No temperature control of room/lab required



FAULT FINDING

Extensive fault finding and testing



STEP-BY-STEP HELP

Screen prompts assist with step-by-step instructions to operate the calorimeter



16 CALIBRATION FIELDS

For different mode and different calorimeters



STANDARD THREAD BOMB VESSEL

Thread bomb vessel



ACCURATE

0.2% RSD accuracy



COMPENSATION*

Compensation for firing energy and sulphur



OPTIONAL BALANCE INTERFACE

Balance interface with baud speed setting



USER FRIENDLY

User Friendly Operation



PRESET FIELDS

One default setting per mode



IMPROVED INITIAL TIMING

Based on drift and time, or time after drift



10 CALIBRATION AVERAGE

For variable amount of calibration average to suit your application



OPERATING PARAMETERS

Operating parameters can be changed via USB interface in experimental mode



FULL LIMS SUPPORT

For the assignment, scheduling, and tracking of samples.



LOW POWER CONSUMPTION

Very low power consumption. No temperature controlling required.



RESTRICT ACCESS

Operating parameter access is password restricted



FILTER DATA ON EVENTS

Extensive testing and detailed data viewing



ECO FRIENDLY

Eco Friendly - small carbon footprint. No water, low power consumption.



LARGE STORAGE

Up to 1024 results storage



TEMPERATURE RANGE

Extensive temperature range from 0°C to 70°C.



INTELLIGENT VESSEL

Intelligent vessel with built-in temperature sensing



EVENT LOGGING

Built-in event logging for ~6000 events



LINEAR SENSORS

Linear temperature sensing with platinum sensors



RESULTS

Results in MJ/Kg, BTU/lb or Cal/g



SAFETY

Safety checks guarantee the safety of the operator.



NO WATER REQUIRED

No Water Bucket. No Spillage. No Measuring.



MANUAL OXYGEN FILLING

Makes use of an external oxygen filling station

CAL3K-S CALORIMETER

The CAL3K Bomb Calorimeter Installation 3K-S-Kit includes :

- Power Supply (External 12V/1.25A)(3K-1-055)
- PC Keyboard PS2 (3K-1-061)
- Preparation Stand (3K-4-49)
- Tweezers (3K-1-081)
- Certified Benzoic Acid Tablets (150 x 0.5g tablets per bottle)(3K-4-71)
- Wire Brush (3K-4-106)
- Printed Installation Manual (INSTALLATION-MANUAL)
- USB 32Gb Green Memory Flash Drive (3K-1-043)

CAL3K-S KIT AND CONSUMABLES

The CAL3K Standard Thread Vessel Installation 3K-5S-Kit includes :


- Complete Top and Bottom Centre Electrode (3K-4-96)
- Outside Electrode (3K-4-37)
- Single Crucibles (3K-4-47)
- Deflector Plate (3K-4-92)
- Firing Cotton (Bundle of 100 threads)(3K-4-65)
- Firing Wire (5 per packet)(3K-4-93)
- Lid O-Ring (3K-4-94)
- Top and Bottom O-Ring in Vessel Lid (3K-4-22)

The CAL3K Filling Station Installation 3K-3-Kit includes :

- Nozzle O-Ring for Defiller Cap and Filling Station (3K-3-29)
- Nozzle O-Ring Used in 3K-3 Jet Assembly (3K-1-080)
- Oxygen Regulator Connection Kit (3K-3-21)
- Defiller Cap (3K-3-22)
- Emergency Deflate Cap (3K-3-18)
- High Pressure Pipe 4mm (Clear/White) (PIPE-4MM-CLEAR)



TECHNICAL SPECIFICATIONS



LCD Display

Large display for easy viewing.



Compact Size



Lightweight

Light weight for easy moving.



TUV CE Certification

Complies with ASTM, DIN and ISO International Standards.

Specification	Information
Working (Operating) Temperature	15-70°C
Storage Temperature	0-70°C
Temperature Resolution	0.000001°C
Reproducibility/Repeatability	0.2% RSD
Resolution	0.0001 MJ/Kg
Results per hour	1 Sample per hour
Measuring range max.	99MJ, 99000J
Working temperature min.	1°C
Working temperature max.	50°C
Temperature Measurement Resolution	10ppm (parts per million)
Cooling Medium	Fan Air

TECHNICAL SPECIFICATIONS

Specification	Information
Type of Cooling	Built-in cooling
Oxygen Operating Pressure Max	40 bar
Balance/Scale Interface	RS232, 1200 to 38400 Baud (settable)
Printer Interface	RS232, 1.2Kb to 115.2Kb
Power Input	2.4W
Interface External Keyboard	PS2
Oxygen Filling	Manual
De gasification	Manual
Halogen (Decomposition) Vessel	Yes, optional
Analysis according to DIN 51900	Yes
Analysis according to ASTM D240	Yes
Analysis according to ASTM D4809	Yes
Analysis according to ASTM D5865	Yes
Analysis according to ASTM E711	Yes
Analysis according to ISO 1928	Yes
Permissible Ambient Temperature	1-35°C
Permissible Relative Humidity	80%
RS232 Interface	Yes
Voltage	220-240 / 100-120V, 12VDC, 1Amp
Frequency	50/60 Hz

Please Note : Technical Specifications subject to change without prior notice.

Please contact our team for accurate technical specifications at the time.



SYSTEM COMPARISON

FEATURE	CAL3K-AP	CAL3K-A	CAL3K-S	CAL3K-F
LIMS	2 x via RS232, or Bluetooth	Yes	Yes	Yes
MASS HEAP	No	No	No	No
BALANCE INTERFACE	Yes, from 1.2 to 38.4KB	Yes, from 1.2 to 38.4KB	Yes, from 1.2 to 38.4KB	Yes, from 1.2 to 38.4KB
RESULT MEMORY	1024 records, 262KB	1024 records	300 records	300 records
TEMPERATURE RESOLUTION	0.000'001°C	0.000'001°C	0.000'001°C	0.000'001°C
DISPLAY	4 x 40 character LCD	4 x 40 character LCD	2 x 20 character LCD	4 x 40 character LCD
KEYBOARD	QWERTY, External, PS2	QWERTY, External, PS2	QWERTY, External, PS2	QWERTY, External, PS2
SAMPLE ID	16 characters, auto-increment	16 characters, auto-increment	16 characters, auto-increment	16 characters, auto-increment
GROUP ID	16 characters	16 characters	-	16 characters
VESSEL RECORD	Yes, unlimited	Yes, unlimited	Limited	Limited
REAL TIME	Yes	Yes	Yes	Yes
CHASSIS IDENTIFICATION	Yes, number	Yes, number	Yes	Yes
CALIBRATION	Normal average & PC	Normal average & PC	Normal average & PC	Normal average & PC
HISTORY CALIBRATION	Yes, up to 10 runs	Yes, up to 10 runs	Yes, up to 16 runs	Yes, up to 10 runs
UNITS	KJ, BTU, CAL	KJ, BTU, CAL	KJ, BTU, CAL	KJ, BTU, CAL
RESULT COMPENSATION	Via PC (IntelCal), Default	Via PC (IntelCal), Default	Yes	Yes
RESULT VALIDATION	Yes	Yes	Yes	Yes
VESSEL PRESS. MONITOR	Up to 100 bar	No	No	No
OXYGEN FILLING	Internal, automatic filling	External manual filling station	External manual filling station	External manual filling station
DE-FILLING	Automatic	Manual	Manual	Manual
MAX CHASSIS RECORDING	Yes	Yes	Yes	Yes
CHASSIS NAME	20 characters, Bluetooth name	20 characters, Bluetooth name	-	-
PASSWORD	CAL3K & PC Password	CAL3K & PC Password	CAL3K & PC Password	CAL3K & PC Password
NETWORK MULTIPLE CALORIMETERS	No	No	No	No
VESSEL LEAK MONITOR	Yes, flags result and warning	No	No	No
EVENT STORAGE	~6000 events	~6000 events	-	~3000 events
EVENT TYPES	~70 different events	~70 different events	-	~80 different events
EVENT CLASSIFICATION	Operational & Technical	Operational & Technical	-	Yes

SYSTEM COMPARISON

FEATURE	CAL3K-AP	CAL3K-A	CAL3K-S	CAL3K-F
VESSEL LOCKOUT, LOCK-IN	Yes, Manual/Auto Linking	Yes, Manual/Auto Linking	Yes, Manual/Auto Linking	Yes, Manual/Auto Linking
SAMPLE REPEAT SPEED (FIRST 2 SAMPLES)	6 min	4-5 min	30 min	7-8 min
OPERATOR TIME PER TEST	+/- 3 min	+/- 3 min	+/- 3 min	+/- 3 min
COOLING	Air	Air	Built-in	Air
COOLING MODES	Ambient/Fixed	Ambient/Fixed	Ambient/Fixed	Ambient/Fixed
RSD	<0.1	0.1	0.2	0.1
POWER CONSUMPTION	0-264 VAC 12W	0-264 VAC 12W	12 VAC 6W	12 VAC 6W
POWER SUPPLY	External 12V	External 12V	External 12V	External 12V
WATER CONSUMPTION	None	None	None	None
REPEATABILITY	<0.1%	0.1%	0.2%	0.1%
CALORIMETER TYPE	Dynamic, Isothermal, Adiabatic	Dynamic, Isothermal, Adiabatic	Dynamic	Dynamic
NUMBER OF VESSELS	Unlimited	Unlimited	Limited (1)	Unlimited
CLOSURE TYPE	Bayonet Lid	Bayonet Lid	Screw Cap (thread)	Bayonet Lid
TESTS P/H WITH 2 VESSELS	10+	10+	1-2 Sample with 1 Vessel Only	8+
BOMB VESSEL TYPE	Removable	Removable	Removable	Removable
OXYGEN FILLING	Automatic	Manual	Manual	Manual
BOMB VESSEL WASHING	Manual	Manual	Manual	Manual
PRINTER CONNECTION	RS232	RS232	RS232	RS232
BALANCE CONNECTION	RS232	RS232	RS232	RS232
ENVIRONMENTAL	5-40°C	5-40°C	5-40°C	5-40°C
PRINTING OF RESULTS	Via PC Software	Via PC Software	Via PC or RS2232	Via PC or RS2232
PC SOFTWARE	Advanced	Advanced	Advanced	Advanced
CORRECTION FACTORS	8	8	4	8
MASS ENTRY	Auto & Manual	Auto & Manual	Auto & Manual	Auto & Manual
CE/TUV CERTIFICATE	Yes (Pending)	Yes (Pending)	Yes (Pending)	Yes (Pending)
VESSEL DETERMINATIONS	Unlimited	Unlimited	5000	Unlimited
SPIKING	Yes	Yes	Yes	Yes
SELF TESTING	Yes	Yes	Yes	Yes

SYSTEM COMPARISON

FEATURE	CAL3K-AP	CAL3K-A	CAL3K-S	CAL3K-F
CONNECTIVITY	USB 2.0, 2 x RS232 at 115.2KB for Bluetooth	USB 2.0, 2 x RS232 at 115.2KB for Bluetooth	USB 2.0, 2 x RS232 at 115.2KB	USB 2.0, 2 x RS232 at 115.2KB
STATS	Yes	Yes	Yes	Yes
PRINTING	Yes, D1 port, 1.2 to 115.2KB	Yes, D1 port, 1.2 to 115.2KB	Yes	Yes
MOISTURE COMPENSATION	Yes	Yes	Yes	Yes
FOOD FIBRE COMPENSATION	Yes	Yes	Yes	Yes
LIMS	Yes	Yes	Yes	Yes
RESULT APPROVAL	Yes, keyboard or PC	Yes, keyboard or PC	Yes, keyboard or PC	Yes, keyboard or PC
REAL TIME PRINTOUT	Yes, optional customer and parameter header	Yes, optional customer and parameter header	-	Yes, optional customer and parameter header
GELATINE CAPSULE	Yes	Yes	Yes	Yes



The CAL3K-S Oxygen Bomb Calorimeter System can be used with most applications including, but not limited to : Animal Feed Research, University Research, Food/Nutrition Analysis, Explosives Analysis, Coal Analysis, Oil Analysis, and other traditional and non-traditional applications.

For more details and application notes visit our website at www.ddscalorimeters.com

CONTACT US

COMPANY HISTORY

Digital Data Systems (DDS has more than 40 years of experience in calorimetry.

In 1972, DDS produced their first calorimeter, the AMPC (Automatic Micro Processor Calorimeter). The AMPC was a dual water isothermal unit controlled by a microprocessor.

In 1980 work began on a new revolutionary design of vessel, namely the DRY vessel or CP510, which meant that there was no surrounding water jacket. A copper sleeve pressed over the vessel replaced the water jacket and the temperature sensors were placed inside the vessel resulting in the heat transfer being extremely fast. Determination time was significantly reduced, increasing the unit efficiency by 4 times. With the processing power of the microprocessors available at the time, the CP500 Calorimeter was born. The striking "buttercup yellow" colour gave a splash of brightness to the then drab laboratories.

In 2002 work began on the CAL2K. The tried and tested DRY system was retained and only the very latest electronic technology was used, including the surface mount devices.

In 2005, DDS came to realize the need for smaller, low volume, inexpensive calorimeter systems, with the same accuracy and reliability of the CAL2K. The ECO was then created as an alternative system to the CAL2K. The ECO is suitable for the following markets: Universities, Research Facilities, Brick Manufacturers, Animal Feed Industries, Food Quality, and Food Production.

In 2007 the new E2K system was developed.

In 2014 work began on the new CAL3K range. The CAL3k-A the first of the range, the CAL3k-F and CAL3K-AP shortly thereafter. In 2018 work began on the new CAL3K-S.

digital data systems (pty) ltd

dds

**DDS Calorimeters are proudly manufactured by :
Digital Data Systems (Pty) Ltd.**

For more information about any of our products visit
our website at www.multifloinstruments.co.in

Contact Us @
Email : sales@multifloinstruments.co.in
Tel : 91-9867349370

Raheja District II Tesla
Unit 142
Navi Mumbai, India
400705

